


## Background

A unique autologous cellular therapeutic (SUPLEXA) has been developed from human PBMC. It is comprised of NK cells,  $\gamma\delta$  T cells and CD8+ T effector cells, capable of broadly lysing a variety of tumor cell lines in vitro. SUPLEXA cells are manufactured using an efficient 2 weeks xeno-free manufacturing procedure employing two proprietary engineered leukocyte stimulator cell lines (ENLIST) that express an array of immunomodulatory proteins. The SUPLEXA cell manufacturing process is highly reproducible and demonstrates low inter-subject variability in cellular composition. SUPLEXA cells are distinguished from many other cellular approaches in that they are derived from autologous PBMC that have only been stimulated with ENLIST cells through naturally occurring receptors without any genetic modification. **Here, we present our *in vitro* manufacturing process for SUPLEXA cells from PBMCs with comprehensive immunophenotyping of SUPLEXA cells.**

The SUPLEXA cells manufacturing process will use peripheral blood mononuclear cells (PBMCs) from cancer patients. PBMCs are stimulated with ENLIST cells for a 5-day induction period, which is then followed by a 9-day cytokine-induced expansion period. SUPLEXA cells are then cryopreserved to use as an autologous adoptive immunotherapy. A first-in-human clinical trial for this adoptive cellular immunotherapy for cancer is projected to this year.



## Methods

**ENLIST cells:** Engineered SK-MEL2 melanoma cell lines (APX-DC and APX-L) that express curated sets of > 20 different immunomodulatory proteins that are engineered for membrane expression. ENLIST cells were used as a lyophilized cellular induction reagent for SUPLEXA.

**SUPLEXA:** PBMCs from normal healthy volunteers were co-incubated with freeze/thaw dead ENLIST cells for 5 days followed by expansion in culture medium containing cytokine support. After 9 days, SUPLEXA cells were harvested and cryopreserved.

**Mass Cytometry (CyTOF):** SUPLEXA cells were comprehensively characterized by mass cytometry (CyTOF) using a 47-marker antibody panel. CyTOF data analysis was done using OMIQ for dimensional reduction by opt-SNE and cell subset phenotyping.


**Tumor Cell Killing Assay:** Tumor cytolytic activity was measured by flow cytometry using fluorescent tumor cell targets at 2:1, 1:1, and 1:2 effector:target cell ratios.

**Cytokines:** A 33 cytokine Luminex panel was used to assess cytokine levels in tumor cell cytotoxicity supernatants.

## Results

**Figure 1: SUPLEXA Cell Manufacturing.** ENLIST immunomodulatory cells are mixed with PBMCs and cultured for 5 days to activate PBMCs. Activated cells are then expanded for 9 days in IL-7 and IL-15. Photomicrographs of 5-day activated PBMCs are shown along with a table listing broad tumor cell killing activity.

### ENLIST Cell Induction of SUPLEXA from PBMCs




### SUPLEXA Cells Develop Broad Tumor Cell Killing Activity


Unmodified Cell Lines Used in Secondary Effector Assay	PBMC coincubated with Unmodified SK-MEL-2 (% killing)	PBMC coincubated with ENLIST Cells* (% killing)
<b>Melanoma Cell Lines</b>		
• SK-MEL-2	0%	>95%
• M14	0%	>95%
• SK-MEL-28	0%	>95%
<b>Non-Melanoma Cell Lines</b>		
• COLO205 - colorectal	0%	>95%
• HOP62 - lung	0%	>95%
• IGR-OV1 - ovarian	0%	>95%
• SKOV - ovarian	0%	>95%
• PC3 - prostate	0%	>95%
• SN12C - renal	0%	>95%
• K562 - leukemia	0%	>95%
<b>Autologous Leukemia Cells</b>	0%	>95%

**Figure 2: Tumor Cell Cytotoxicity Assay and Killing Activity.** SUPLEXA cells were compared for tumor cell cytotoxic activity against fluorescent tumor target cells using a flow cytometry method. Figures show representative killing of M14-RFP target cells as measured by flow cytometry and tumor cytotoxicity activity against different tumor cell lines (M14, K562, PC3, COLO205).

### Flow Cytometry Tumor Killing Assay; SUPLEXA:M14-RFP Targets




### Comparison of SUPLEXA Tumor Cell Line Killing Activity



**Figure 3: Cytokine Production by SUPLEXA During Tumor Cytotoxicity.** SUPLEXA cells were incubated without or with M14, K562, PC3, or COLO205 cells for 48 hours. Supernatants were tested for cytokines by 33-plex Luminex. Plots of tumor-induced cytokine production by SUPLEXA are shown.

### Cytokine Production by SUPLEXA Corresponds to Cytotoxicity Activity




## Results


**Figure 4: SUPLEXA Cell Phenotyping By CyTOF – Antibody Panel and Data Analysis Workflow.** SUPLEXA cells were analyzed to compare single-cell phenotypes using a customized 48-marker CyTOF antibody panel to identify NK cell and T cell subsets as well as transcription factors, adhesion markers, and functional molecules. CyTOF antibody panels, CyTOF data analysis workflow, and a cell subset overlay plot is shown here.

Marker	Antibody Clone	Label
CD45	HI31	Pan-leukocyte marker
CD172ab	SE8A5	Mesoderm cell marker
CD8a	RPA-78	CD8 T cell alpha-chain
CD20	2H7	B cell marker
CD4	RPA-74	CD4 T cell marker
CD3	UCHL1	AB T cell marker
CD11b	DX11	Adhesion molecule for NK cell activation
Granzyme A	Polyclonal	Effector cell tumor killing
KIR2DS1	1127B	NK receptor that interacts with DAP12
PD-1	EHF23H7	Activator and checkpoint inhibitor
DAP12	406288	Signaling adaptor of NK cell activation
CD19	HB19	B cell marker
CD14	HEB2	monocyte and macrophage
CD86	IT2-2	costimulatory
Granzyme B	CB9	Effector cell tumor killing
Granzyme M	GM26E7	Effector cell tumor killing
CD20	B4D3	Dispositional neutrophil-associated disease
TCR- $\gamma$	REA173	gamma delta T cell receptor
CD134	LD10	Checkpoint inhibitor
CD29	PNB3	Lymphocyte activation marker
TCR- $\delta$	REA591	gamma delta T cell identification
CD28	WD1028	Transcription factor for effector lymphocyte function
CD8	SD388E	CD8 T cell beta-chain
CD16	304	NK marker
CD39	A1	ATP to Adenosine enzyme
PC2	R17-439	Transcription factor, regulates NK cells
PC3	P39-19	NK cell marker
PSGL1	CHD131	Adhesion molecule for tissue infiltration
CD56	NCA196.2	Pan-NK cell marker
4-1BB	REA108	T cell costimulation
GITR	110416	Glucocorticoid-induced receptor - Tregs
NKG2A	131411	Effector NK cell marker - negative
CD137a	HA63	Coaggregation marker
Grz B	GB11	Effector cell tumor killing
NK4F8	REA808	Activating receptor on NK cells
TCR- $\beta$	REA179	Invariant NK T cells
CD137	2H6-2A3	PD-1 ligand for checkpoint inhibition
TCR- $\gamma$	12383	gamma delta T cell receptor
NKG2D	148910	Effector NK cell marker - positive
TCR- $\alpha$	T1881-18-31	alpha beta T cell receptor
Ki67	8D5	Recently proliferated cell marker
Thal1	4B10	Transcription factor for NK and effector lymphocyte subsets
Helios	4929	Effector cell tumor killing
CD57	REA789	T cell exhaustion marker
HLA-DR	L243	Incompatibility
SLAMF6	202811	Trigger cytotoxic activity of NK and NKT cells
CD38	HT2	NK adhesion to endothelium, cAMP diase hydrolyase
TIGIT	A15153G	Inhibitory receptor on Tregs

### CyTOF Data Analysis Workflow




### Cell Subset Overlay opt-SNE Plot




**Figure 5: SUPLEXA Cell Phenotyping By CyTOF.** CyTOF data from SUPLEXA were analyzed to compare functional phenotypes by marker expression profiles.

### Cell Identification Marker Expression



### Functional Marker Expression



## Conclusions

1. The manufacture of SUPLEXA cells from PBMCs generates immune effector cells with potent cytotoxicity activity against multiple tumor cell line targets.
2. Comprehensive phenotyping of SUPLEXA cells by CyTOF indicates that SUPLEXA cells are a mixture of NK cells, CD8+ T cells, CD56+ NK-like T cells, and  $\gamma\delta$  T cells that express activation phenotypes (NKG2D, CD26) and high-levels of killer cell factors (granzymes and perforin).
3. SUPLEXA cells are a novel autologous cellular therapeutic for cancer.